


Annexure 1: Percolation Vessel (Capacity - 100LTRS)




Percolation Vessel (Capacity - 50LTRS)



Detail of Label Pocket






- 1.Tight lid
- 2.vessel need to be fixed in the round clamp, vessel should be removable and rotateable one.
- 3.bottom should be provided with SS ball valve with nipple at the end for tube connetion
- 4.Latch removable and lockable setup required to place a vessel at the bottom
- 5. Lockable trolley wheel need to be provided
- 6. SS Round mesh with 5 MM mesh hole
- 7. SS bracket with food grade seal for close packing

	ifecare Limited CORPORATE HIL BHAVAN INMENT OF INDIA ENTERPRISE HILL SHAVAN KERALA - 69:	E & REGD. OFFICE N, POOJAPPURA NTHAPURAM	Percolation Vessel 50 Ltrs	
DRAWN BY:	NK		LILL /CLIO /DDO L/EVT/LIOMCO /OOOE O/	REVISION: 00
CHECKED BY:	I R			DATE:18.09.25
APPROVED BY:	harikrishnan k p		SITE: HOMCO New Factory, Velya Kalavoor	ISSUED FOR:
DATE:	12.09.2025			APPROVAL
SCALE: ALL DIMENSIONS	NOT TO SCALE ARE IN MILLIMETRES		This drawing is the Sole Property of M/s.HLL Lifecare Ltd. The contents of this drawing shall not be reproduced in any manner or used for other purposes the Mind intended one.	SHEET NO:

Annexure - III Collection vessel 50 liters



Annexure - IV Collection vessel 100 liters

